

### **MBO GUI User Guide Overview**





Multi-ROI, multi-Zplane statistics

### **Quickly Process with Suite2p**





### Save Full Dataset to .tiff, .zarr, .bin, or .h5



## Installation

- Create environment UV
- Install mbo\_utilities >= 2.0.1
   (Packaged with cellpose, suite2p, suite3d, graphical utilities)
- 3. Start main GUI: uv run mbo






# **Data Selection**

- General data loader, opens supported filetypes (.tiff, .zarr, .bin, .h5)
- "Data Preview" widget is available for Raw ScanImage Tiffs
- Non ScanImage tiffs will open as a 3D or 4D array


Here, we focus on the *Data Preview widget* 


Hover (?) for tooltips, many buttons also contain tooltips on hover



# **Open File / Select Folder**

**Open Files:** Dialog box that shows selectable and multi-selectable





# **Open File / Select Folder**


**Select Folder**: Dialog box that allows folder selection only

- Will not show you files in that folder
- Helpful for large datasets that may contain multiple files

There are files here, as seen previously.

In this case, all files in this directory will be loaded.

Opens tiff / zarr / bin / h5 produced by mbo.imwrite()




# **Load Options**

Only affect data visualization (e.g. no mROI splitting happens on disk)

### (Raw ScanImage tiffs only) Separate ScanImage mROI's


- Open ScanImage multi-ROI's as separate arrays
  - Visualizations (e.g. scan-phase correction) is done PER-ROI
  - Helpful to compare scan-phase offsets between mROI's
- If unchecked, both mROI's share an optimal phase-offset value that is calculated when frames are vertically concatenated (as they are on disk, pre-stitching)
- If checked, mROI's have distinct phase-offset values

# Unchecked (stitch/concatenate mROIS)





# Checked (Separate mROIS)





### **Load Options**

Only affect data visualization (e.g. no mROI splitting happens on disk)

### **Metadata Preview (experimental)**

- Quickly preview all available consolidated and raw scanimage metadata
- Supports raw ScanImage tiffs, and any extracted outputs from mbo.imwrite
- Also support suite2p ops.npy files





### Collapsable subgroups


```
si.version 4
    ▶ si.FrameData
    ▼ si.RoiGroups
       ▼ si.RoiGroups.imagingRoiGroup
         si.RoiGroups.imagingRoiGroup.ver 1
si.RoiGroups.imagingRoiGroup.classname 'scanimage.mroi.RoiGroup'
         si.RoiGroups.imagingRoiGroup.name 'MROI Imaging ROI Group'
          si.RoiGroups.imagingRoiGroup.roiUuid '6594A3E141FEAC46'
          si.RoiGroups.imagingRoiGroup.roiUuiduint64 7.319655482e+18
           ▼ si.RoiGroups.imagingRoiGroup.rois [list]
              ▼ si.RoiGroups.imagingRoiGroup.rois[0
                si.RoiGroups.imagingRoiGroup.rois[0.ver 1
                si.RoiGroups.imagingRoiGroup.rois[0.classname 'scanimage.mroi.Roi'
                si.RoiGroups.imagingRoiGroup.rois(0.roiUuiduint64 1.218622949e+19 si.RoiGroups.imagingRoiGroup.rois(0.zs [0]
                  ▼ si.RoiGroups.imagingRoiGroup.rois[0.scanfields
                   si.RoiGroups.imagingRoiGroup.rois[0.scanfields.ver 1
                   si.RolGroups.imagingRolGroup.rois[0.scanfields.roiUuid '9353A1919E060303' si.RolGroups.imagingRolGroup.rois[0.scanfields.roiUuiduint64 1.061600639e+19
                   si.RoiGroups.imagingRoiGroup.rois[0.scanfields.centerXY [-3.814104639, 0]
                   si.RoiGroups.imagingRoiGroup.rois[0.scanfields.sizeXY [7.213114754, 18.03278689]
                    si.RoiGroups.imagingRoiGroup.rois[0.scanfields.enable 1
                    ► si.RoiGroups.imagingRoiGroup.rois[0.scanfields.pixelToRefTransform [list]
                    ➤ si.RoiGroups.imagingRoiGroup.rois[0.scanfields.affine [list]
                si.RoiGroups.imagingRoiGroup.rois[0.discretePlaneMode 0
                si.RoiGroups.imagingRoiGroup.rois[0.enable 1
        ► si.RoiGroups.integrationRoiGroup
size 2453342976
```

## **Preview Data Widget (Raw ScanImage Tiff Only)**

- Lazy-load a raw-scanimage tiff
- Larger tiffs (e.g. 60,000x14x400x400) could take ~30s to load on the first load
- After that, it will be cached and load in sub 5s


#### **Window Functions**

- View mean, max, std, mean-subtracted image over the given window size
- Windows > 20 frames will slow rendering
- Windows are applied over the time dimension only
- Windows do not take effect until Window Size > 3 frames



### sigma

- Apply a spatial gaussian filter on each frame
- Helps visualize very common pre-processing step supplied as a parameter to processing pipelines
  - suite2p "smooth sigma" parameter for registration




## **Preview Data Widget (Raw ScanImage Tiff Only)**

- Lazy-load a raw-scanimage tiff
- Larger tiffs (e.g. 60,000x14x400x400) could take ~30s to load on the first load
- After that, it will be cached and load in sub 5s


#### **Window Functions**

- View mean, max, std, mean-subtracted image over the given window size
- Windows > 20 frames will slow rendering
- Windows are applied over the time dimension only
- Windows do not take effect until Window Size > 3 frames

- Dynamic range can vary drastically between projections
- Use the 'auto-adjust contrast' button to automatically adjust vmin/vmax



alculate vmin/vmax using the full dataset (can be computationally intensive for large datasets)



You can also adjust the vmin/vmax manually with the histogram widget

## **Preview Data Widget (Raw ScanImage Tiff Only)**


- Lazy-load a raw-scanimage tiff
- Larger tiffs (e.g. 60,000x14x400x400) could take ~30s to load on the first load
- After that, it will be cached and load in sub 5s

#### **Window Functions**

- View mean, max, std, mean-subtracted image over the given window size
- Windows > 20 frames will slow rendering
- Windows are applied over the time dimension only
- Windows do not take effect until Window Size > 3 frames
- Dynamic range can vary drastically between projections
- Use the 'auto-adjust contrast' button to automatically adjust vmin/vmax



alculate vmin/vmax using the full dataset (can be computationally intensive for large datasets)



You can also adjust the vmin/vmax manually with the histogram widget

(Raw ScanImage Tiff Only)

#### **Scan Phase Correction**

- Compare with / without phase correction
- Use a projection with window-size of 3-20
  - Large datasets can be very slow with window size> 6
- Scan-offset results are typically reliable with a window size > 3

Fix-Phase: Enable / Disable scan-phase correction

Sub-Pixel: Enable / Disable FFT-based sub-pixel phase correction

Upsample (sub-pixel only): lines are registered within 1/upsample of a pixel

- Larger upsample = more precise shifts = more computationally intensive


Exclude border px: Pixels on the border tend to be more poorly aligned due to the non-linear angular velocity of the resonant galvo scanner. This excludes a number of pixels around the full image or individual mROIs.

*Max-offset:* Limit the allowed offset to this many pixels. Helpful in noisy recordings when some frames with minimal signal yield offsets 5+. Offsets very rarely exceed 2 or 3 pixels.

### Recommended workflow:

View a mean-subtracted 3-15 window-size projection.

Check all mROI's without Fix Phase checked, and with fix-phase checked. First adjust border/max offset, see if you can get visually reliable correction. Turn on and off sub-pixel rapidly (may need to lower window size) and see if it improves scan-phase correction. If so, decrease upsample until you find a stable result across many frames (typically a value of 2 or 3).



(Raw ScanImage Tiff Only)

### **Summary Stats**

- Compute mean/std/snr for each mROI regardless of splitting or stitching on every 10<sup>th</sup> frame
- SNR: mean / std\_img
- Combined: Take mean of all mROI's

Useful for evaluating:


- z-plane quality as a function of depth
- mROI quality over the FOV

\*note: these means are used for meansubtracted image preview

▶ N N □ ✓ S framerate: 20

▶ K N □ ✓ S framerate: 20

- + fps




(Raw ScanImage Tiff Only)

### **Saving Data to Disk**

- File -> Save As
- Opens a transparent save dialog

\*Note: This demo is fast because it's a single ~1GB tiff

I/O speed depends on file-size, how users chose to split frames across tiffs (in scanimage during acquisition), output filetype, FFT and upsampling factor.




(Raw ScanImage Tiff Only)

#### **Saving Data to Disk**

- File -> Save As
- Opens a transparent save dialog

\*Note: This demo is fast because it's a single ~1GB tiff


I/O speed depends on file-size, how users chose to split frames across tiffs (in scanimage during acquisition), output filetype, FFT and upsampling factor.



(Raw ScanImage Tiff Only)

### **Reading Extracted Data**

- Ensure you TURN OFF data preview or this will error
- This will not yet allow split-mROI merging



### **Suite2p Processing**

- Run Suite2p on *currently selected z-plane*
- All parameters, with descriptions, provided
- Use crop selector to run on a spatially cropped subset of the data

